

# Pytentiostat documentation

Oliver Rodriguez

## Module: **potentiostat**

Module to select the potentiostat to use and perform electrochemical experiments (CV, LSV, CA, OCP).

### **pytentiostats.Setup(model=0, path\_exe='', folder='', verbose=1)**

Initializes the potentiostat to use and selects folder to save data files.

#### **Parameters**

- **model:** *str* or *False*. Potentiostat model to use. Accepted values: 0, False, 'chi760e'
- **path\_exe:** *str*. Directory where the potentiostat executable is located.
- **folder:** *str*. Name of an existing directory where data files will be saved.
- **verbose:** *boolean*. Print information on the console. 1 by default.

#### **Methods**

*pytentiostats.Setup.info()* Prints the potentiostat model, executable path and directory to save the data files.

### **pytentiostats.Technique(text='', fileName='CV')**

Class to initialize electrochemical techniques.

#### **Parameters**

- **text:** *str*. String containing the macro or commands to send to the potentiostat.
- **fileName:** *str*. Base filename to save data files.

#### **Attributes**

- **bpot:** *boolean*. False by default. If True, activates the second working electrode.

#### **Methods**

- *Technique.writeToFile()*. Saves the macro or commands sent to the potentiostat as a file.
- *Technique.run()*. Performs the measurement and calls writeToFile to save the macro.
- *Technique.message()*. Used to output to console when the experiment starts and finishes.
- *Technique.bipot(E=-0.5, sens=1e-6)*. If *pytentiostats.Technique.bpot* is set to True, it will activate the second working electrode with the specified parameters.

### **pytentiostats.CV(Eini=-0.2, Ev1=0.2, Ev2=-0.2, Efin=-0.2, sr=0.1, dE=0.001, nSweeps=2, sens=1e-6, fileName='CV', header='CV')**

#### **Parameters**

- **Eini:** *double*. [V] initial potential.
- **Ev1:** *double*. [V] first vertex potential.
- **Ev2:** *double*. [V] second vertex potential.
- **Efin:** *double*. [V] final potential.
- **sr:** *double*. [V/s] scan rate.
- **dE:** *double*. [V] potential increment.
- **nSweeps:** *int* number of sweeps.

- sens: *double*. [A/V] current sensitivity.
- fileName: *str*. Name to save data file. File will be overwritten if it exists.
- header: *str*. Header of the txt data file. Use it to include comments.

**pytentiosats.LSV(Eini=-0.2, Efin=0.2, sr=0.1, dE=0.001, sens=1e-6, fileName='LSV', header='LSV')**

#### Parameters

- Eini: *double*. [V] initial potential.
- Efin: *double*. [V] final potential.
- sr: *double*. [V/s] scan rate.
- dE: *double*. [V] potential increment.
- sens: *double*. [V/A] current sensitivity.
- fileName: *str*. Name to save data file. File will be overwritten if it exists.
- header: *str*. Header of the txt data file. Use it to include comments.

**pytensiostats.CA(Estep=0.2, dt=0.001, ttot=2, sens=1e-6, fileName='CA', header='CA')**

#### Parameters

- Estep: *double*. [V] potential step to apply.
- dt: *double*. [s] time increment.
- ttot: *double*. [s] total time of the step.
- sens: *double*. [V/A] current sensitivity.
- fileName: *str*. Name to save data file. File will be overwritten if it exists.
- header: *str*. Header of the txt data file. Use it to include comments.

**pytentiosats.OCP(ttot=2, dt=0.01, fileName='OCP', header='OCP')**

#### Parameters

- ttot: *double*. [s] total time of the step.
- dt: *double*. [s] time increment.
- fileName: *str*. Name to save data file. File will be overwritten if it exists.
- header: *str*. Header of the txt data file. Use it to include comments.

## Module: file

Performs common operations on files created by the pytensiostat library.

**pytentiosats.Read(text=False, model=False)**

Reads a file created with the pytensiostat library.

#### Methods

**\*pytentiosats.Read.read(text=0, model=0)**

- text: *str* or *False*. Text to locate in file to use as starting point to read.
- model: *str* or *False*. Potentiostat model used when creating the files.

**\*pytentiosats.Read.search(text)**

- text: *str*. Text to locate in file to use as starting point to read.

**pytentiostats.LoadXY(fileName='file', folder='.', skiprows=0, delimiter=',', model=0)**

Loads a general data file with X and Y columns.

**Parameters**

- *fileName*: *str*. Filename of the file to read with extension included.
- *folder*: *str*. Directory where the file is located.
- *skiprows*: *int*. Number of rows from the top to skip.
- *delimiter*: *char*. Character used to delimiter columns: ',', '\t', etc.
- *model*: *str* or False. Potentiostat model used when creating the files.

**Returns** data: LoadXY instance - x: *numpy array* - y: *numpy array*

**pytentiostats.LoadCV(fileName='file', folder='.', model=0)**

Loads a cyclic voltammogram data file obtained with the pytentiostat library.

**Parameters**

- *fileName*: *str*. Filename of the file to read with extension included.
- *folder*: *str*. Directory where the file is located.
- *model*: *str* or False. Potentiostat model used when creating the files.

**Returns** data: LoadCV instance

- *E*: *numpy array*. [V] potential array.
- *i*: *numpy array*. [A] current array.

**pytentiostats.LoadLSV(fileName='file', folder='.', model=0)**

Loads a linear sweep voltammogram data file obtained with the pytentiostat library.

**Parameters**

- *fileName*: *str*. Filename of the file to read with extension included.
- *folder*: *str*. Directory where the file is located.
- *model*: *str* or False. Potentiostat model used when creating the files.

**Returns** data: LoadCV instance

- *E*: *numpy array*. [V] potential array.
- *i*: *numpy array*. [i] current array.

**pytentiostats.LoadCA(fileName='file', folder='.', model=0)**

Loads a chronoamperogram data file obtained with the pytentiostat library.

**Parameters**

- *fileName*: *str*. Filename of the file to read with extension included.
- *folder*: *str*. Directory where the file is located.
- *model*: *str* or False. Potentiostat model used when creating the files.

**Returns** data: LoadCV instance

- *t*: *numpy array*. [s] time array.
- *i*: *numpy array*. [i] current array.

**pytentiostats.LoadOCP(fileName='file', folder='.', model=0)**

Loads a chronoamperogram data file obtained with the pytentiostat library.

**Parameters**

- *fileName*: *str*. Filename of the file to read with extension included.

- folder: *str*. Directory where the file is located.
- model: *str* or False. Potentiostat model used when creating the files.

**Returns** data: LoadCV instance

- t: *numpy array*. [s] time array.
- E: *numpy array*. [V] potential array.